

# Capacitance Point level detection Nivector FTC968

For all types of powdered and fine-grained solids



## Lợi ích:

- Simple mounting and commissioning without calibration
- Mechanical safety, cost-efficiency and long operating life due to no wearing parts
- Reliable operation due to high immunity to electromagnetic interference and voltage peaks
- Simple control thanks to the visibility of the switching status from outside the vessel
- Function test also with filled silo possible

từ **111,00 €**

Giá kể từ 23.04.2021

**Phạm vi ứng dụng:** Nivector FTC968 is a simple and cost-effective probe for point level detection in bulk solids in silos containing free-flowing powders or fine-grained solids. Its flush-mounted construction and materials make it suitable for use in the food industry.

Tín năng và thông số kỹ thuật