

Nivector FTC968



Lợi ích:

- Simple mounting and commissioning without calibration
- Mechanical safety, cost-efficiency and long operating life due to no wearing parts
- Reliable operation due to high immunity to electromagnetic interference and voltage peaks
- Simple control thanks to the visibility of the switching status from outside the vessel
- Function test also with filled silo possible

Phạm vi ứng dụng: Nivector FTC968 is a simple and cost-effective probe for point level detection in bulk solids in silos containing free-flowing powders or fine-grained solids. Its flush-mounted construction and materials make it suitable for use in the food industry.

Tín năng và thông số kỹ thuật