

Memosens Analog Converter CYM17

Converter for Memosens sensors for easy use with laboratory fermenters



Thông tin thêm và mức tính giá hiện tại:

www.apsc.endress.com/CYM17

Lợi ích:

- Easy exchange of your existing analog sensors. Just install the Memosens sensors to your fermenter and connect them with the converter.
- Various adapter cables are available for easy connection of one Memosens pH sensor (e.g. CPS171D) and one optical DO sensor (e.g. COS81D).
- Thanks to the Memosens technology, the cable connection is insensitive to moisture, corrosion or dirt: you can count on the sensor's highest availability and your trial is safe.
- Memosens sensors have an excellent performance in process applications as well as in benchtop fermenters. Providing you with 100% measuring consistency from the first lab trials to the final scaled-up process.
- The converter continuously checks the Memosens sensors for highest measuring reliability and sets off an alarm if any error occurs.

Phạm vi ứng dụng: The Memosens Analog Converter CYM17 enables the easy use of digital Memosens sensors in your fermentation applications in the laboratory without changing your SOP. Memosens digital technology combines maximum process and data integrity with simple operation. It eliminates any moisture issues and guarantees stable signals after autoclavation.

Tín năng và thông số kỹ thuật

pH

Measuring principle

Potentiometric

pH

Application

CYM17 converter enables easy retrofit of conventional analog pH and oxygen sensor to Memosens CPS171D and COS81D in benchtop fermentation applications.

Material

Housing: Aluminum

Dimension

Height: 36.6 mm (1.44 in)

Deep: 72 mm (2.8 in)

Width 140.6 mm (5.5 in)

Process temperature

-5 to 50 °C (20 to 120 °F)

Ingres protection

IP54

Input

M12 socket

Output / communication

Output pH: mV signal according to Nernst equation with 0 mV at pH 7 and -59 mV/pH

Temperature pH: simulated PT1000

Output oxygen: 0 to 60 nA with 60 nA at 100% oxygen saturation

Temperature oxygen: simulated 22 kOhm NTC

Oxygen

Measuring principle

Amperometric oxygen measurement

Oxygen

Application

CYM17 converter enables easy retrofit of conventional analog pH and oxygen sensor to Memosens CPS171D and COS81D in benchtop fermentation applications.

Material

Housing: Aluminum

Dimension

Height: 36.6 mm (1.44 in)

Deep: 72 mm (2.8 in)

Width 140.6 mm (5.5 in)

Process temperature

-5 to 50 °C (20 to 120 °F)

Input

M12 socket

Output / communication

Output pH: mV signal according to Nernst equation with 0 mV at pH 7 and -59 mV/pH

Temperature pH: simulated PT1000

Output oxygen: 0 to 60 nA with 60 nA at 100% oxygen saturation

Temperature oxygen: simulated 22 kOhm NTC

Thông tin bổ sung www.apsc.endress.com/CYM17